




Silicone Tube Loop

May 2021

A silicone loop can be added to any 8610C SRI GC as long as there is a place for the 10port valve which is required. A multiple Gas #5 GC configuration is shown in the photo.

In this GC there are three 10 port valves in the heated valve oven.

The front panel oven of the valve oven is easily disconnected by removing six 5/32" hex head screws.



SRI Tech Support: 310-214-5092

www.srigc.com

SiliconeTubeLoopMay2021

Page 1




Silicone Tube Loop

May 2021

Ideally the 3rd valve will be located at the front where the loop connections (ports 3 and 10) are easiest to get to.

Feed the silicone loop through the front plate and attach to ports 3 and 10 of the valve. If there is a regular loop already attached to those ports you will need to remove it in order to attach the silicone tube loop.

The silicone loop consists of two 1/16" stainless steel tubes with a short section of silicone tubing attached to the ends. The stainless tubing is .010 diameter to minimize the volume of the stainless tubing. .010" tubing has a volume of about 1.3ul/inch so two 24" long stainless tubes have a total volume of about 64 ul



SRI Tech Support: 310-214-5092

www.srigc.com

SiliconeTubeLoopMay2021

Page 2



Silicone Tube Loop

May 2021

The ends of the stainless tubing are epoxied to a plastic barb fitting. These fittings are commonly available in various plastics or metals. Silicone tubing is attached to the barb fittings.

The tube shown here is 12 inches long and has an internal diameter of .058 which has a volume of 43ul/ inch for a total of 516 ul. Longer or shorter tubing can be substituted as desired.

The silicone part of the loop is inserted into the liquid to be measured.



SRI Tech Support: 310-214-5092

www.srigc.com

SiliconeTubeLoopMay2021

Page 3

Silicone Tube Loop



May 2021

The silicone tubing can be any size or wall thickness, but we have used this which was available on-line.

The valve diagram shows how Valve3 either isolates the loop (no gas flow) or places the loop in the flow path of valve2.

As the silicone loop sits in the liquid many of the compounds in the liquid will permeate through the silicone wall into the interior of the loop. When valve2 rotates and then Valve3 also rotates, whatever is in the silicone loop is pushed into the column.

Valve 3 is rotated back to the Load position (loop isolated) after a few seconds to avoid a detector signal shift which occurs if the carrier gas continues to extract the loop after it is initially emptied.



SRI Tech Support:
310-214-5092


SiliconeTubeLoopMay2021

Page 4

Silicone Tube Loop

May 2021

As an example, the silicone tube loop (MultipleGAs#5 GC with TCD and FIDmethanizer) was alternately immersed in 10% methanol/water mix and clean water. This is the chromatogram of the methanol. You can see some small CO2 peak in the beginning also.



These are the operating conditions.

The most important point is that the analysis has to be a repetitive sequence exposing the silicone tube to the sample for a constant amount of time each run.

In this case the runtime was 3.9 minutes and the Postrun Delay before re-starting was .1 minutes. So there was a 4 minute cycle.

The results of many sequential runs were automatically added to the log file which is shown on the next page.

Channel 1 events
E:\SiliconeLoop\SiliconeLoop.evt

Time	Event
0.000	ZERO
0.020	F ON (Valve2Rotate)
0.040	E ON (Valve3Rotate)
0.300	E OFF (Valve3Rotate)
1.750	INTEG IMMEDIATE
2.900	INTEG IMMEDIATE
3.000	F OFF (Valve2Rotate)

Channel 1 temperature control
E:\SiliconeLoop\Silicone.tem

Init temp	Hold	Ramp	Final temp
170.00	3.900	0.000	170.00

178.50
0.00
3.900

Channel 1 post-run actions

- Save file as: SiliconeFID11.CHR Auto-increment
 - Or use list of filenames: List
- Save results Use data file name Use fixed file name: _____
- Add to results log: CH1.LOG
- Print results Update DDE link Save picture
- Execute: _____
- Restart run after: 0.10 minutes 0 times total (0 remaining)
- Recalibrate at level: 1 Save results file to FTP site
- Smooth first Save data file to FTP site
- Copy data to channel: 2
- Add to 3D display
- Generate signature
- Match signatures
- Email:
 - No
 - On alarm condition
 - Always



SRI Tech Support: 310-214-5092

www.srigc.com

SiliconeTubeLoopMay2021

Page 5

Silicone Tube Loop


May 2021

Here is the log file which shows the data from each run as the sample was switched from 10% Methanol/Water to Clean Water and then back again to 10% Methanol/Water.

Note the Relative Standard Deviation (%RSD) for the first 7 runs was 1.27%.

Also note how the silicone tube takes several runs to equilibrate to the new sample. This may take longer with more viscous fluids or other target molecules.

For this reason, you have to sample on some repeated interval so the sample has the same amount of time every time to permeate through the wall of the silicone tubing.



```
*CHI - Notepad
File Edit Format View Help

Start test by placing silicone loop in 10% methanol/water
andriy88.CHR 5/4/2021 12:26:55 "Methanol-FID" 2.160 3065.7094
andriy89.CHR 5/4/2021 12:30:55 "Methanol-FID" 2.166 3089.2204
andriy90.CHR 5/4/2021 12:34:56 "Methanol-FID" 2.170 3027.6563
andriy91.CHR 5/4/2021 12:38:56 "Methanol-FID" 2.173 2983.7288
andriy92.CHR 5/4/2021 12:42:56 "Methanol-FID" 2.173 3009.3998
andriy93.CHR 5/4/2021 12:46:56 "Methanol-FID" 2.160 3048.5872
andriy94.CHR 5/4/2021 12:50:56 "Methanol-FID" 2.163 3031.9966

Calculate %Relative Standard Deviation for above 7 runs =1.27%

Switch to clean water. Change data file names
SiliconeFID01.chr 5/4/2021 12:54:56 "Methanol-FID" 2.083 249.7427
SiliconeFID02.chr 5/4/2021 12:58:56 "Methanol-FID" 2.086 14.9428
SiliconeFID03.chr 5/4/2021 13:02:56 "Methanol-FID" 2.100 3.4962
SiliconeFID04.chr 5/4/2021 13:06:56 "Methanol-FID" 0.000 0.0000

Notice how it takes several runs to reach zero.

Switch back to 10% methanol in water
SiliconeFID05.chr 5/4/2021 13:10:56 "Methanol-FID" 2.090 1.1803
SiliconeFID06.chr 5/4/2021 13:14:57 "Methanol-FID" 2.203 2772.3704
SiliconeFID07.chr 5/4/2021 13:18:57 "Methanol-FID" 2.183 3092.3008
SiliconeFID08.chr 5/4/2021 13:22:57 "Methanol-FID" 2.183 3067.3421
SiliconeFID09.chr 5/4/2021 13:26:57 "Methanol-FID" 2.180 3048.2049
SiliconeFID10.chr 5/4/2021 13:30:57 "Methanol-FID" 2.183 3125.6104
SiliconeFID11.chr 5/4/2021 13:34:57 "Methanol-FID" 2.180 3055.5135
```



SRI Tech Support: 310-214-5092

www.srigc.com

SiliconeTubeLoopMay2021

Page 6